



A NEW GENERATION OF THE iRT KIT

Our bestseller updated.

NEW iRT KIT GENERATION LAUNCHED: SAME KIT, BETTER PERFORMANCE

Biognosys' iRT Kit is a peptide kit specifically designed to be spiked in every analyzed sample for LC-MS quality control and retention time normalization. After many years of successful application, we have updated our bestseller to support you even better in these tasks.

SAME PEPTIDES, SAME PHILOSOPHY, SAME WORKFLOWS

The same 11 stable, non-naturally-occurring peptides elute at the same times in your LC setup. You won't need to adjust any workflow involving our kit.

BALANCED SIGNALS

A more balanced intensity of peptide signals minimizes the issue of undetected peptides under suboptimal LC conditions.

QUANTITATIVE ACCURACY

The iRT Kit peptides are now AAA quantified, fully enabling their use for quantitative quality control.

Are you interested in the QC capabilities of the iRT Kit?
Unlock its full potential by using it with our **QuiC software**.

For further details visit:



[www.biognosys.com/shop /iRT-kit
/QuiC](http://www.biognosys.com/shop/iRT-kit/QuiC)